# MEDIAWEL 10

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## English version

➔ Use of the machine .....	48
➔ Use of the documentation .....	48
■ Training obligation.....	48
■ Available documentation .....	48
■ Use of the operating manual.....	48
➔ Assistance.....	49
■ Troubleshooting .....	49
■ Spare parts .....	49
➔ Warranty .....	49
<b>Safety and the environment .....</b>	<b>50</b>
➔ General safety instructions.....	50
➔ Risks and prevention .....	50
➔ Environment .....	51
■ Discharging waste water .....	51
■ Scrapping the machine.....	51
<b>Installation - Commissioning .....</b>	<b>52</b>
➔ Unpacking.....	52
■ Unpacking procedure.....	52
■ List of supplied accessories .....	52
➔ Environment requirements.....	53
■ Overall dimensions.....	53
■ Site requirements .....	53
➔ Connection .....	54
■ Water.....	54
■ Power supply.....	54
➔ Initial commissioning .....	54
■ Switching on .....	54
<b>Presentation .....</b>	<b>55</b>
➔ Description of the device .....	55
■ Front .....	55
■ Lid details .....	55
■ Bucket inside.....	56
■ Rear .....	56
➔ Presentation of the control screen .....	57
■ Home screen .....	57
■ Lights.....	57
➔ Screens navigation .....	58
➔ Operating principle .....	59
<b>Parameters .....</b>	<b>60</b>
➔ Settings.....	60
■ Date and time .....	60

■	Language .....	60
➔	<b>Operators .....</b>	<b>61</b>
■	Screen.....	61
■	Create.....	62
■	Modify and delete.....	62
➔	<b>Programs.....</b>	<b>63</b>
■	Screen.....	63
■	Create.....	64
■	Modify and delete.....	70
➔	<b>Maintenance.....</b>	<b>70</b>
	<b>Manage the access rights.....</b>	<b>71</b>
➔	<b>Users.....</b>	<b>71</b>
➔	<b>Activation.....</b>	<b>71</b>
➔	<b>Audit trail.....</b>	<b>71</b>
	<b>Cycle.....</b>	<b>72</b>
➔	<b>Use .....</b>	<b>72</b>
➔	<b>Cycle screen.....</b>	<b>72</b>
➔	<b>Cycle operating .....</b>	<b>74</b>
■	Preparation.....	74
■	Programming and starting .....	75
■	Delayed start.....	77
■	Interruption .....	78
■	Alarms and icons .....	78
■	Distribution.....	80
■	Cleaning .....	81
➔	<b>Stirring .....</b>	<b>82</b>
■	Principle .....	82
■	Screen.....	82
➔	<b>E-trace traceability .....</b>	<b>83</b>
■	USB .....	83
■	Online service Labpage .....	83
➔	<b>Servicing and maintenance.....</b>	<b>84</b>
■	Summary table .....	84
■	Temperature probe calibration.....	84

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ALLIANCE BIO EXPERTISE, Guipry-Messac, France

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
The screenshot shows the control interface for the MEDIAWEL 10 system. At the top, the logo for ALLIANCE BIO EXPERTISE is on the left, and the model name 'MEDIAWEL 10' is on the right. The main display is divided into two columns. The left column, titled 'CYCLE', lists menu options: PROGRAMME, AGITATEUR, EAU, PARAMETRAGE, and DROITS D'ACCES. Below these is a large digital display showing '20.1 °C'. At the bottom of the left column are icons for user profile, system status, and a refresh button. The right column displays the current cycle: '01: TEST'. It lists two steps: '1: STERILISATION' with parameters '121.0°C - 015' - 9 - 9' and 'No Alarm', and '2: DISTRIBUTION' with parameters '42.0°C - 000' - 9 - 3' and 'No Alarm'. At the bottom of the interface are navigation buttons: left arrow, right arrow, a green checkmark button, and a grey 'X' button. The bottom left corner shows the date '03/10/2019', time '10:56:52', and version 'V2.0.01/A008'.

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




The screenshot shows the 'MEDI AWEL 10' interface. On the left, a menu lists 'CYCLE', 'PROGRAMME', 'AGITATEUR', 'EAU', 'PARAMETRAGE', and 'DROITS D'ACCES'. Below the menu, the temperature is set to '16.0 °C'. At the bottom left, there are icons for user, printer, and other functions, along with the date '18/12/2019' and time '17:31:30'. The main area displays '01: TS' and 'Démarrage différé' with date and time selection fields set to '18/12/2019' and '17:31'. A green checkmark and a red 'X' icon are visible at the bottom right.

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The screenshot shows the 'MEDI AWEL 10' interface. On the left, the menu is the same as in the previous screenshot. The temperature is now '16.2 °C'. The main area displays '01: TS' and a graph showing a temperature profile. The graph is labeled 'EN PAUSE' and shows a temperature of '0.0' and '16.2 °C'. At the bottom, a yellow bar contains the text 'Démarrage du cycle à' followed by a black box with '22/12/2019 17:31' and a red 'CANCEL' button. The date and time at the bottom left are now '18/12/2019 17:34:02'.

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
The image displays two screenshots of the MEDIAWEL 10 interface. Both screens show the 'ALLIANCE BIO EXPERTISE' logo and the 'MEDIAWEL 10' title. The left screenshot shows the login screen with a 'Password' field and a numeric keypad. The right screenshot shows the maintenance menu with 'TextLabel' fields for 'Real temperatures (reference probe)' and 'Measured temperatures (internal probe)', each with 'Low step' and 'High step' input fields. A green checkmark button is visible in the bottom right corner of the right screenshot.

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**ENGLISH VERSION**

**Declaration of Conformity MEDIAWEL30**
**ALLIANCE BIO EXPERTISE, Guipry-Messac, France**

Declares that the here below material complies with the following european directives:

Description	Models
<b>MEDIAWEL 10</b>	<b>MEDW1010</b>

Comply with

**EU directives:**

Pressure Equipment	2014/68/EU
Electromagnetic Compatibility	EN61326-1/2013
Safety requirements for electrical equipment for laboratory use	IEC61010-1: 2010
	EN610010-2-10:2014

Material resistance calculation according :

**ASME VIII div 1 Edition 2013**
**ASME II Part D Edition 2013**
**Directive 2014/68/EU**

<b>Description of pressure vessel :</b>	
Fluid group	2
Category (2014/68/EU)	1
Max operating pressure	1,4 bar
Safety valve set pressure, Max allowable pressure PS	1,7 bar
Allowable temperature TS	0-123 °C
Volume	18 l
Test pressure (water)	1,7 bars
Dessin No./Rev	140405-21-100rev05

Declaration issued by the manufacturer

The person responsible of the declaration is:

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<b>Made in FRANCE</b>			
<b>ALLIANCE BIO EXPERTISE</b>			
230~10% V~	3600 W	16Amp.50-60 Hz	Weight: 63Kg
<b>Pressure vessel identification: N°</b>		190014-3	
Drawing number: 180405-1-1 rev 4			
Fluid group	2	Allowed temperature TS	0-123°C
Category (2014/68/EU)	1		
Max operating pressure	1,2 bar	Test pressure PT	1,72 bar
Safety valve setting & PS	1,7 bar	Pressure test media	water
<b>Serial N°: 519.....</b>		<b>N° : MEDW1010</b>	

## ➔ Use of the machine

The MEDIAWEL must only be used if it is in perfect condition and in strict compliance with the instructions set out in this manual.

Safety glasses must be worn when using the device.

## ➔ Use of the documentation

### ■ Training obligation

Use of this kind of equipment by untrained and inexperienced personnel may pose a risk to the user and people in its vicinity.

Users must be trained to operate the MEDIAWEL by an authorized distributor or the manufacturer.

Users must have been made aware of the risks associated with the operation of pressurized steam devices.

### ■ Available documentation

A maintenance manual is available and placed at the disposal of experienced technicians, provided they have completed technical training provided by the manufacturer.

### ■ Use of the operating manual

The operating manual must have been read in its entirety before using the MEDIAWEL for the first time. It must be consulted should an error message be displayed.

The preventive maintenance operations set out in this manual must be carried out at the frequency stipulated in this manual.

Account must be taken of the instructions displayed on screen, the pictograms and other materials, otherwise there is a risk of injury to users and damage to equipment.

## Recipients

The following tasks are authorised :

- **User**
  - Creation of preparation cycle programs,
  - Preparation of the culture media and starting cycles,
  - Stopping cycles and dispensing culture media,
  - Cleaning,
  - Level 1 preventive maintenance operations set out in the manual.
  
- **Maintenance Engineer**
  - Calibration of the temperature probe,
  - Level 1 preventive maintenance operations set out in the manual.

## ➔ Assistance

### ■ Troubleshooting

In the event of a problem, please contact your authorized distributor. State the model and the product's serial number and describe the problem.

We strongly recommend that you take out a maintenance contract for optimum use of the product and for your safety. For all questions relating to maintenance contracts, please contact your authorized distributor or the manufacturer.

ALLIANCE BIO EXPERTISE  
ZA de Courbouton-Le Tremplin  
35480 GUIPRY-MESSAC France

Tel: +33 2 40 51 36 10

e-mail: [service@abioexpertise.com](mailto:service@abioexpertise.com)

### ■ Spare parts

For all spare parts orders, please contact your authorized distributor or the manufacturer.

ALLIANCE BIO EXPERTISE  
ZA de Courbouton-Le Tremplin  
35480 GUIPRY-MESSAC France

Tél : +33 2 40 51 36 10

Email : [service@abioexpertise.com](mailto:service@abioexpertise.com)

## ➔ Warranty

ALLIANCE BIO EXPERTISE hereby certifies that the device displays the characteristics set out in this manual on delivery. Machines are guaranteed against manufacturing defects and defective components for a period of twelve months from invoicing date.

ALLIANCE BIO EXPERTISE shall carry out possible repairs of components, within its premises, during the warranty period. All components or devices experiencing any kind of malfunction must be returned to the ALLIANCE BIO EXPERTISE factory at the cost of ALLIANCE BIO EXPERTISE's representative or the client, if necessary.

The warranty shall not be valid for issues relating to transport, accidents, incorrect use or internal interference caused by a poor quality electricity network or people who are not authorized to use the device.

## Safety and the environment

### ➔ General safety instructions

	<p><b>WARNING</b></p> <p>The MEDIAWEL is exposed to high temperatures and pressures during sterilization cycles. For this reason, it is essential to follow the safety instructions in this operating manual.</p>
	<p><b>STOP</b></p> <p>In case of emergency, unplug the wall power supply before any manipulation.</p> 

- Never use the automated preparator to sterilize hazardous materials or products, which contain alkalis. Sterilization of these products may cause an explosion, corrosion of the chamber or chamber pipework and damage to seals.
- Should a liquid containing saltwater and high salinity agar-saline agar, etc. be used, you must rinse the vessel well after use and wipe the latter's surfaces as well as around the cover seal (corrosion risk).
- Do not modify the device.
- The MEDIAWEL must form the subject of annual preventive maintenance, carried out by trained personnel authorized by ALLIANCE BIO EXPERTISE.
- Defective parts must only be replaced with original spare parts supplied directly or indirectly by ALLIANCE BIO EXPERTISE.

### ➔ Risks and prevention

Symbol	Risk	Cause	Component	Means of prevention	Instructions
 	Burns	Hot surface	Lid Bowl	Thermal insulation of components Lid and hatch locked during the working cycle No access above 80°C and if pressure >0	Do not touch hot surfaces Wear gloves
		Splashing	Dispensing and adding port caps	No access above 80°C and if pressure >0	Release the steam before opening the lid by using the safety valve
	Burns	Water steam	Safety valve Caps	No access above 80°C	Wear gloves
	Pinching	Lid weight	Lid	Counter balance springs	Hold the lid during closing

## ➔ Environment

### ■ Discharging waste water

Water drained from the cooling system can be discharged into the rainwater system.

The system must be designed to accommodate the discharge of boiling water.

Water used to clean the vessel must be discharged into the wastewater system.

### ■ Scrapping the machine

This device contains metal, plastic, electrical and electronic components.

Directive 2002/96/EC on the recovery and recycling of Waste Electrical and Electronic Equipment (WEEE) requires Electrical and Electronic Equipment (EEE) to be recycled at the end of its life.

For this reason, end users must dispose of electrical and electronic waste at dedicated waste centers or collection points.



## Installation - Commissioning

### ➔ Unpacking

#### ■ Unpacking procedure

This chapter covers unpacking the device and lists the accessories supplied.

- Cut the pallet strapping
- Lift the cardboard sleeve and remove it
- Remove the upper foams
- Remove the plastic protection
- By pair lift it by the sides and put it on its place of use.


#### ■ List of supplied accessories

Part No.	Quantity	Name
MEDW2100	1	User manual
MEDW4104	1	Magnetic stirrer
MEDW2103	1	Sampling tube
MEDW4012	4 m	Inlet and outlet water tubing + 4 clamping clips
MEDW3003	1	Bactystoppers bag (x50)
MEDW4100	1	Stirrer bearing kit

## → Environment requirements

### ■ Overall dimensions

### ■ Site requirements

#### Atmosphere

Operating temperature : from 15 to 30°C

Hygrometry: maximum 80% relative humidity, non-condensing

#### Water

Pression : 2 bars minimum

Temperature : 5°C to 20°C

#### Power supply

Tension : 230 V +/- 10 %

Intensity : 16 A

Frequency : 50/60 Hz

Power consumption : 3600 W +/- 10 %


## ➔ Connection

It is essential to fit the power cable with a plug that is suitable for the above power.

### ■ Water

#### Inlet

- For more convenient use, include an easy to access tap close to the MEDI AWEL, at a height of approximately 1 meter and slightly offset at the side of the device as illustrated below
- Cut the 4 meter hose to the desired lengths for supply and outlet hoses
- Fit the in-line filter (ref: MEDW3008) to the water supply hose
- Connect the supply and outlet hoses as set out below

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#### Outlet

The cooling water drainage system must be designed to withstand a temperature of 130°C. It may be made from copper or PVC-C or any other material that is suitable for this use.

### ■ Power supply

The adapted plug must be used adapted to the above specifications :

- **Phase** wire is identified by **1**, printed on one of the 2 black cables
- **Neutral** wire is identified by **2**, printed on the second black cable
- **Earth** wire is yellow and green

The wall power supply must be between 1 and 1,2m height and equipped with a shutdown system.

## ➔ Initial commissioning


### ■ Switching on

Prior to switching the device on, you must ensure that the mains voltage between phase and neutral complies with 230V +/- 10%. Plug the machine using the plug into the wall power supply and switch on the front switch.

# Presentation


## ➔ Description of the device

### ■ Front


1	Lid
2	Handle for opening the lid
3	USB Port
4	Master switch
5	Safety cover
6	Touch screen
7	Vessel pressure gauge
8	Draining tray location


### ■ Lid details

1	Vessel safety valve
2	Sampling port
3	Adding port
4	Cover seal
5	Handle for unlocking the lid
6	Cotton wool housing


■ Bucket inside

1	Bucket
2	Temperature probe
3	Stirrer

■ Rear


1	Supply hose connection
2	Drain hose connection
3	RJ45 plug (Ethernet)

## ➔ Presentation of the control screen

### ■ Home screen

The screen is organised as follows:


The upper section of the left menu is used for navigation in the menus in standby mode, and to display the temperature and the remaining time for a step, in cycle mode.

The lower section of the left-hand column enables the date and time to be displayed. The main display area allows ongoing actions to be displayed.

The « validate » and « cancel » keys are displayed in the lower right-hand corner of the screen. Other dynamic keys are displayed depending on specific requirements.


### ■ Lights

There is a light on each side of the control panel. Its role is to indicate, alert or inform the user of the status of the current cycle, of an alarm or a status. It indicates this by means of a light that is on or off, flashing or permanently on, and is white, orange, red, blue or green, based on the rules appearing in the following table:

Light color	Status	Description of the status
White	Permanently on	Stand-by
	Flashing	
Blue	Permanently on	Processing to Media preparation
	Flashing	
Orange	Permanently on	Warning message during cycle (Safety cover or doortrap not closed)
	Flashing	Warning message during cycle (step message, temperature out of tolerance)
Red	Permanently on	
	Flashing	Major default during cycle (heating default, cooling default)
Green	Permanently on	Media ready for dispensing
	Flashing	Waiting for the acknowledgement of the safety message "Make sure that vessel pressure is 0 bar before opening"

## ➔ Screens navigation

## ➔ Operating principle

The automated preparator is a device for the automatic preparation of culture media ("media"), which are used in microbiology laboratories for conducting analyses.

It comprises a vessel equipped with a heating device (electrical heating elements around the vessel), a cooling device (circulation of cold fluid around the vessel), a stirring device inside the vessel and a temperature-measuring device within the medium.


The medium, contained in the vessel, is made up of water and soluble substances.

The automated preparator is equipped with independent safety features (safety valve, thermal circuit breaker, etc.) and safety features controlled by the automated system:

A set of sensors wired in series, enabling the closing of all covers to be controlled.


An electrically actuated bolt enabling the covers to be locked in a closed position

When the temperature measured within the medium is equal to or greater than the *safety temperature*, the covers prohibiting access to hot surfaces must be closed.

# Parameters

## ➔ Settings

### ■ Date and time


When the screen is displayed, the fields are already filled with current timestamping values. Pressing on a text field activates the numeric keypad for this field.

### ■ Language

The arrows scroll through the available languages for the text field. Press on the menu “Cycle” will allow the full upgrade of the screen in the selected language.

## → Operators

### ■ Screen


### Active keys:

CYCLES	Switches back to the first screen in the program choice sequence
PARAMETERS	Inactive
PROGRAMS	Switches to the program parameters screen
USERS	Inactive
MAINTENANCE	Switches to the "maintenance" screen
Arrows	Switch from one operator position to another
Validate	For available positions : access to creation of an operator For already created operators : access to modification and deletion of an operator
Cancel	Not activated

## ■ Create

Using the arrows, chose the next available position and press Validate. Le keyboard appears and write down the operator's name. The yellow arrow erase the last character, the red cross erase the whole name and the green arrow validates the name. Validate, the operator's name is saved.

The screenshot shows the 'MEDI AWEL 10' interface. On the left, there is a menu with options: CYCLE, **PARAMETERS**, PROGRAMS, USERS, SETTINGS, MAINTENANCE, and ACCESS RIGHTS. Below the menu, the temperature is displayed as **15.4 °C**. At the bottom left, there are icons for a person, a gear, a globe, and a document, along with the date and time: 23/12/2019 12:08:43 and version V2.0.0 / A010.

The main area is titled 'Operator N°' and contains a text input field. Below the input field is a virtual keyboard with the following layout:


1	2	3	4	5	6	7	8	9	0	←
A	Z	E	R	T	Y	U	I	O	P	×
Q	S	D	F	G	H	J	K	L	M	↵
W	X	C	V	B	N		.	-	/	↵

## ■ Modify and delete

Please refer to the Modify or Deletion of programs.

## ➔ Programs

### ■ Screen

Active keys:

CYCLES	Switches back to the first screen in the program choice sequence
PARAMETERS	Inactive
PROGRAMS	Inactive
USERS	Switches to the user parameters screen
SETTINGS	Switches to the date and time settings screen
MAINTENANCE	Switches to the "maintenance" screen
Arrows	Scroll through the programs
Enter	For blank programs: switches directly to the "program creation/modification" sequence For existing programs: switches to a choice between the program creation / modification or deletion sequence
Escape	Switches back to the first screen in the program choice sequence

Note : for created programs, the parameters sum up is displayed in the main zone. For empty programs, the program name is replaced by « Available », and the main zone is empty.

## ■ Create

### • Principle

The automated preparator enables a "preparation cycle" to be automatically performed for a volume of culture medium. A preparation cycle comprises  $n$  "steps";  $n = 2, 3$  or  $4$

Each step is defined by a temperature setpoint (in °C) and a duration (in minutes)


The temperature setpoint for a step is equal to or greater than 25.0°C and equal to or lower than 123.0°C. It is defined to the nearest 0.1°C

The duration of a step is equal to or greater than 1 minute and equal to or less than 999 minutes. It is defined to the nearest 1 minute.

The duration of the last step is always infinite (the cycle is stopped by the user).

The start of the cycle is triggered either by the user (manual start), or by means of a timer (scheduled start).

The phases between the steps are transition phases, which may be "heating" or "cooling"

Each program contains the following information:

Field title	Number
Program number	01 to 50
Program name	40 characters
Number of steps $n$	2, 3 or 4
Name of step no. 1	40 characters
Temperature setpoint for step no. 1	25.0 to 123.0°C
Duration of step no. 1	1 to 999
Alarm at the beginning of step no. 1	yes/no
Alarm message for step no. 1	40 characters (if an alarm is enabled)
Stirring speed in transition to step no. 1	1 to 9
Stirring speed during step no. 1	1 to 9
Name of step no. 2	40 characters
Temperature setpoint for step no. 2	25.0 to 123.0°C
Duration of step no. 2	1 to 999 or infinite if $n=2$
Alarm at the beginning of step no. 2	yes/no
Alarm message for step no. 2	40 characters (if an alarm is enabled)
Stirring speed in transition to step no. 2	1 to 9
Stirring speed during step no. 2	1 to 9



The choice of program number is determined by the operation prior to starting the "program creation / modification" sequence.

When creating a program, the fields are empty, so valid values must be entered in all the fields.

When modifying a program, the fields are already populated with the program's current parameters. The user chooses which fields to modify.

- **Creation**

Using the arrows, chose the next « available position » and press «Validate. The screen « Program name » appears, press on the text field to make the keyboard appear. The yellow arrow is to delete a character, the red cross is to delete the whole name, and the green arrow is to validate the name. Press on Validate to save the name and then press on Validate to go to the next screen. Pressa gain on the text fields if necessary.

ALLIANCE  
BIO EXPERTISE

# MEDIAWEL 10

CYCLE

**PARAMETERS**

PROGRAMS

USERS

SETTINGS

MAINTENANCE

ACCESS RIGHTS

---

15.6 °C

18/12/2019  
16:43:36

V2.0.0 / A010

04: TSA

**Number of steps (2 - 4)**

2

ALLIANCE  
BIO EXPERTISE

# MEDIAWEL 10

CYCLE

**PARAMETERS**

PROGRAMS

USERS

SETTINGS

MAINTENANCE

ACCESS RIGHTS

---

15.6 °C

18/12/2019  
16:45:48

V2.0.0 / A010

04: TSA

**Step name 1**

STERILISATION

ALLIANCE  
BIO EXPERTISE

# MEDIAWEL 10

CYCLE

**PARAMETERS**

PROGRAMS

USERS

SETTINGS

MAINTENANCE

ACCESS RIGHTS

---

15.6 °C

18/12/2019  
16:47:41

V2.0.0 / A010

04: TSA

**Step temperature 1 (25.0 - 125.0)**

25.0

ALLIANCE  
BIO EXPERTISE

MEDIAWEL 10

CYCLE

**PARAMETERS**

PROGRAMS

USERS

SETTINGS

MAINTENANCE

ACCESS RIGHTS

---

15.6 °C

18/12/2019  
16:48:33

V2.0.0 / A010

04: TSA

Step duration 1 (1 - 999)

1

ALLIANCE  
BIO EXPERTISE

MEDIAWEL 10

CYCLE

**PARAMETERS**

PROGRAMS

USERS

SETTINGS

MAINTENANCE

ACCESS RIGHTS

---

15.6 °C

18/12/2019  
16:49:13

V2.0.0 / A010

04: TSA

Step alarm 1

Message

NB: The "enabled" changes status each time it is pressed. It contains "enabled" or "disabled" information. The selected key appears in green. If the alarm is disabled, the "message" field is greyed out and is not editable. It is only editable when the alarm is enabled.

ALLIANCE  
BIO EXPERTISE

MEDIAWEL 10

CYCLE

**PARAMETERS**

PROGRAMS

USERS

SETTINGS

MAINTENANCE

ACCESS RIGHTS

---

15.6 °C

18/12/2019  
16:52:57

V2.0.0 / A010


04: TSA

Stirring speed before step 1 (1 - 9)

7

ALLIANCE BIO EXPERTISE MEDIAWEL 10

CYCLE  
**PARAMETERS**  
 PROGRAMS  
 USERS  
 SETTINGS  
 MAINTENANCE  
 ACCESS RIGHTS

---

**15.7 °C**

18/12/2019  
 16:54:39

V2.0.0 / A010

04: TSA


Stirring speed during step 1 (1 - 9)

7

The parameters screens for steps 2 to 4 are identical to those for step 1.

Once all the parameters have been entered, you go back to the list of programs.

ALLIANCE BIO EXPERTISE MEDIAWEL 10

CYCLE  
**PARAMETERS**  
 PROGRAMS  
 USERS  
 SETTINGS  
 MAINTENANCE  
 ACCESS RIGHTS

---

**15.6 °C**

18/12/2019  
 16:52:57

V2.0.0 / A010

04: TSA

Stirring speed before step 1 (1 - 9)

7

- **How to set up the stirrer speed ?**

Standards media

- 2 to 10L

Speed 7 (for heating/ Sterilization/ Cooling)

Speed 4 (for dispensing)

- 1 to 2L

Speed 5 (For heating / Sterilization / Cooling)

Speed 3 (for dispensing)

Viscous media

- 2 to 10L

Speed 7 (for heating / Sterilization/ Cooling)

Speed 5 (for dispensing)

- 1 to 2L

Speed 6 (for heating / Sterilization / Cooling))

Speed 4 (for dispensing)

Low thermal property media

- 2 to 10L

Speed 7 (for heating / Sterilization/ Cooling)

Speed 7 (for dispensing)


- 1 to 2L

Speed 6 (for heating / Sterilization / cooling)


Speed 6 (for dispensing)

## ■ Modify and delete

When an existing program is selected, press on Validate and a menu enabling a choice to be made between modification and deletion of the program appears:

If you want to modify the program, press on “Modify the program”, the Program Name screen will appear. If you want to delete the program, press on « Delete the program”, the following confirmation screen will appear.

Press on “Confirm” to cancel the deletion, or press “Delete the program” to delete the program.


## ➔ Maintenance

This menu is dedicated to servicing and maintenance operations, to be realized by the technical team of Alliance Bio Expertise and/or trained by the manufacturer. Please refer to Servicing and maintenance part.

## Manage the access rights

### → Users

This function is under development.

ALLIANCE BIO EXPERTISE MEDIAWEL 10

CYCLE  
PARAMETERS  
**ACCESS RIGHTS**  
USERS  
ACTIVATION  
AUDIT TRAIL

01: DPAC

16.0 °C

24/12/2019 10:21:27

V2.0.0 / A010

### → Activation

This function is under development.

ALLIANCE BIO EXPERTISE MEDIAWEL 10

CYCLE  
PARAMETERS  
**ACCESS RIGHTS**  
USERS  
ACTIVATION  
AUDIT TRAIL

16.0 °C

24/12/2019 10:22:18

V2.0.0 / A010

Active access rights

Operator	Administrator
<input type="checkbox"/>	<input type="checkbox"/>

### → Audit trail

This function is under development.

# Cycle

## → Use

Users must have familiarized themselves with the previous paragraphs, relating to the environment and safety, before proceeding.

Follow the security procedures :


- Never attempt to force the lid opening.
- Always ensure that the pressure shown by the vessel pressure gauge is 0 before attempting to open the lid.
- Always manually operate the vessel's safety valve before attempting to open the lid, or the adding or decanting port.
- Should you need to move the device between 2 workstations, it is important to move it slowly and to always hold it in both hands during the move.
- Always put on protective gloves designed for high temperatures before accessing the vessel → certain surfaces can actually still be very hot even though the temperature of the medium is below 80 °C.
- Never attempt to force the drain access hatch open.

## → Cycle screen

### • STAND-BY cycle

This screen appears when the device is started in stand-by mode. The last program used is suggested when the device is switched on. Its no., name and a summary of the program's parameters are displayed.

The temperature of the vessel probe is displayed in real time in the left-hand column.


Active keys:

STIRRER		Switches to the "Stirrer" screen
WATER		Force the water circulation in the double jacket (filling or force cooling)-Available in 2018 only.
PARAMETERS		Switches to the "Parameters" screen
Arrows		Increase / decrease the program number and display the corresponding program name
Enter		Confirms the choice of program and moves to the next screen
Escape		Inactive

NB: the box around the program name means that it is a choice from a list.

- **WORKING cycle**





Main area :

- Program number and name
- Display area for the standard curve. The standard curve appears from the beginning of the cycle. The curve background is coloured based on the actual progress of the cycle (on the basis of the transition temperature, on the basis of the remaining time for a step)
- Current phase (based on the program cycle)
- XXX:
  - In transition: setpoint temperature of the next step + "° C"
  - In a step: Remaining time in an HH:MM:SS format
- Medium temperature "° C"
- Alarm zone: displays messages relating to current alarms or errors, with a red background for alarms and errors.

## ➔ Cycle operating


### ■ Preparation

Sequence of operations	Opération	Symbol
1	Connect the water supply hose, Connect the drain hose, Open the water tap	
2	Connect the power cable, enable the system's circuit breaker,	  
3	Make sure the drain plug is securely in place	
4	Fill the hot-water bath with 2.5L of water (Water level just under the 2 <sup>nd</sup> turn of the heating component),	
5	Put the bucket into the vessel, a bucket lock orientated at 10 Go the bucket down to the bottom and turn in the clockwise direction until it is locked in the vessel. <b>PAY ATTENTION NOT TO HIT THE TEMPERATURE PROBE</b>	 Clip
6	Place the stirrer on the axis of the bucket after après s'être assuré de la présence de la rondelle de glisse sur l'axe,	
7	Replace the carded cotton in the cage (Before each cycle)	
8	Switch on the device by pressing the main switch	
9	Pour 1/3 of demineralized water of the total volume to be prepared in the vessel	
10	Slowly pour the powder, avoiding putting it on the wall and on the blades of the stirrer	
11	Start the agitator at speed 5	


12	Pour the additional demineralized water over the wall contour and the stirrer the powder to remove the powder deposit on the metal parts, to reach the total volume to be prepared (1 to 10 liters, 12 maximum if the Mediawel is not moving (no cart)).	
10	Close the lid, and lock it by pushing the handle Close the safety cover and make sure that the LEDs turn blue	
11	Select the program using the navigation arrows and follow the instructions given by the various display that will scroll after pressing the green “validation” button to the start screen (see next paragraph “running a program” for more details).	
12	Press the green button when the message “Cycle start?” appears.	

■ Programming and starting


Once the above steps done, select the program required thanks to the arrows, and press on Validate :


Select the operator using the arrow and press Validate :


Press on the grey field to open the keyboard and write the batch number : press on Validate. Optional field.

Press on the grey field to open the keyboard and write the total volume present inside the bucket (between 1 and 10 liters) and validate.


Start the cycle : check the different security element before pressing on Validate.


### ■ Delayed start

On the start cycle screen, press on the clock, the following screen appears.

Press on the arrows to select the date and time the cycle should start then press validate. The following screen shows the delayed start information, and the possibility to abort it by pressing on Cancel.


## ■ Interruption

The cycle can be interrupted at any stage. Simply pressing the key stops the cycle.

Interrupting a cycle immediately results in cold water flowing, if the temperature of the medium was less than 80 °C.

In all cases, the vessel should only be **accessed after having enabled the vessel's safety valve (pull the ring upwards)** and when wearing protective gloves designed for high temperatures. In effect, surfaces may still be at very high temperatures and can cause burns.

**ALLIANCE BIO EXPERTISE** **MEDIAWEL 10**

**CYCLE**

- PROGRAM
- STIRRER
- WATER
- PARAMETERS
- ACCESS RIGHTS

01: TS

TRANSITION  
STERILISATION  
121.0  
**15.4 °C**

**Cycle stop : please confirm**

OK ANNULER


23/12/2019  
12:12:49

V2.0.0 / A010

## ■ Alarms and icons

### • Alarms

In case of color status change on the machine, a message appears in the « message zone » to inform users about the action to do.

**CYCLE**

- PROGRAM
- STIRRER
- WATER
- PARAMETERS
- ACCESS RIGHTS

## : Program name

Current step XXX  
55,1 °C

**Zone Message**

24/03/2016  
16:46:35

Message	Action
Step Start Alarm	Follow the instruction of the message that have been created for this message. Example: "Add blood"
Alarm 'Close Vessel Lid'	The lid is not properly closed. Release the pressure from the vessel by using the safety while the temperature is below 80°C. Then open the lid and close it again to get blue light coming one. If not →Call service.
Alarm 'Close Safety Covers'	One of the safety lids (door trap or lid to access to the dispensing port) is not closed or not properly closed--> Close the lids
Heating error 'Abnormal temperature variation'	Time to reach the sterilization set point is too long or the temperature is not stable --> Make sure that the stirrer is running properly. If not, check the bearing state.
Cooling error 'Abnormal temperature variation'	Time to reach the dispensing set point is too long or the temperature is not stable --> Make sure that the water tap is open.
Temperature warning 'temperature is out of tolerance'	The temperature was not in +/-1°C range during the sterilization plateau-Make sure that temperature is within the tolerance of your application before validating the cycle.
Restart cycle after power on	Just for information - Check the data records at the end of the cycle before validating the cycle.
'Media' temperature probe error 'Abnormal temperature measurement'	The vessel temperature probe is defective --> Call service
Technical Error 'Stirrer motor'	Defective stirrer --> Call service
Technical Error 'Pressure release solenoid'	Defective pressure release valve --> Call service
Technical Error 'Cooling'	Defective cooling valve --> Call service
Technical Error 'Locker magnet'	Defective lid locking solenoid --> Call service
Communication Initialisation failed	Switch off and on again the machine - Call service if it does not solve the problem
Communication failed. Please switch off and contact the after-sales service	Switch off and on again the machine - Call service if it does not solve the problem

### • Icônes

- Communication between HMI and automaton cards normal

- Failure of communication between HMI and automaton cards normal

- SD card connected and compliant


- SD card failure



- Connected to an Ethernet network




- USB key connected

■ Distribution

Operations order	Operations	Symbols
1	<p>The end of the cycle is indicated by the lights changing from blue to green</p> <p>This means that the dispensing temperature has been reached and at its step → dispensing can be undertaken</p>	
2	<p>Open the safety cover</p>	
3	<p>Pull the ring on the safety valve in order to ensure that any possible pressure is released</p>	
4	<p>Prepare the decanting tubing and place it on top of the lid but leave it in its protective packaging before removing the decanting cap</p>	
5	<p>Heat the adding port with a blow torch before inserting the tip of the decanting tubing</p> <p>Do not heat the decanting tip to not damage the seal</p>	
6	<p>Connect the hose to the pump and run dispensing immediately</p>	

■ Cleaning

Sequence of operations	Operations	Symbols
1	Remove the dispensing hose and immediately rinse it with very hot water	
2	Open the lid and take off the bucket turning it counterclockwise. <b>PAY ATTENTION NOT TO HIT THE TEMPERATURE PROBE</b>	
3	Wash the bucket and the stirrer with very hot water using a soft brush Wash the sampling and addition plugs with very hot water	
4	Place a container under the device	
3	Remove the drain plug until the container is full Drain the entire volume and rinse the vessel with very hot water	
4	Disable the installation circuit breaker Disconnect the plug from the power cord	
5	Remove the carded cotton from its cage and clean it with 70° alcohol	
6	Remove the cover seal and rinse it with very hot water	
7	Leave the lid open until the next use and only fill the hot-water bath before using the appliance Do not put the bucket and drain plug back into the vessel until the next use	

## ➔ Stirring

### ■ Principle


The stirrer's speed can be set between 4 rpm and 64 rpm by selecting a figure between 1 and 9 in the stirring speed parameters.

The following table displays the corresponding values in rpm:

Speed parameter	1	2	3	4	5	6	7	8	9
Corresponding speed in rpm	4 rpm	15 rpm	25 rpm	40 rpm	50 rpm	54 rpm	58 rpm	62 rpm	64 Rpm

The speed setting can be set in stand-by mode and must be configured in each program for all stages of the cycle (see § stirrer in stand-by mode and § program creation) .

### ■ Screen

Active keys:

PROGRAMS	Switches back to the first screen in the program choice sequence without stopping the stirrer
PARAMETERS	Switches to the "Parameters" screen
Arrows	Increase / decrease the stirrer's speed (0 = STOP; 9 = max.) - Immediate application
Enter	Go back to the program selection screen without stopping the stirrer.
Escape	Switches back to the first screen in the program choice sequence without stopping the stirrer


The stirrer speed can be set up from 0 to 9.

## ➔ E-trace traceability

### ■ USB

The Etrace module on serie allow the data export on a USB key. The key has to be connected on the Mediawel **before starting the cycle**. The key can be disconnected 30 seconds after the final validation of the cycle, under .CSV format or PDF files.

The PDF file is a report of the cycle which allow a quick validation of the cycle. The CSV file allow to handle raw information for complementary details about events or status of the device. This file can be required by the after-sales service.

ID Prep	S/N	Batch n	Date/Time	User	Status	IHMV	AUTV									
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:22	AS1	START	v0.3.9	0									
MEDIAWEL30	517016000	170513 1022	13/05/2017 11:22	AS1	END	v0.3.9	0									
ID Prep	S/N	Batch n	Date/Time	Number of step	Temp Step1	Step1 Durati	Step1 Stir sp	Temp Step2	Step2 Durati	Step2 Stir sp	Volume	NPROG	PRNAM			
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:22	2	121.0	15	3	45.0	0	3	3.0		5 AS1 STIRRER			
ID Prep	S/N	Batch n	Date/Time	Status	Media Temp	Default code										
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:22		30	57.3										
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:23		30	57.2										
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:23		30	57.4										
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:23		30	57.0										
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:23	Transition step1		55.6										
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:23	Transition step1		54.6										
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:23	Transition step1		54.2										
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:24	Transition step1		53.9										
MEDIAWEL30	517016000	170513 1022	13/05/2017 10:24	Transition step1		53.6										

### ■ Online service Labpage

The Mediawel is delivered with a 3 months free licence. The Etrace can be connected by WIFI (integrated into the Mediawel), or by Ethernet cable, to access the Labpage service - Real Time Data Monitoring - Historic directly from your PC and smartphone :

- Real Time Data Monitoring : follow in real time the temperature chart, addition and distribution step alarms without being in the lab !
- Historic : up to 10 000 cycles stocked on the Etrace to be downloaded by CSV and PDF for local archive. Can be sorted out on different values (date, user, volume ...) to edit statistics and numeric archives of all the quality information of your media !

**Ask you engineer for the installation !**

## ➔ Servicing and maintenance

### ■ Summary table

Element	Operation	Fréquency				
		1 cycle	20 cycles	2 months	6 months	year
Bactystopper MEDW3003	Replacement	•				
Sampling tube O ring (X5) MEDW4101	Replacement			• Regular check required especially if flamed		
O ring for addition cap MEDW4103	Replacement					•
Lid seal MEDW4102	Replacement					•
Water inlet filter (Option) MEDW3008	Cleaning				•	
Stirrer bearing kit MEDW4100	Replacement					•
Ring seal on stirrer axis MEDW4008	Replacement				•	

### ■ Temperature probe calibration


the manufacturer strongly recommends calibrating the temperature probe once a year.

- Fill the vessel of 10 liters of water
- Close the lid
- Remove the sampling cap from the lid
- Insert reference temperature probe (42 cm long probe with specific adapter)
- Launch a cycle with a program at 121°C of 15 min speed 4 for the first step and 42°C speed 7 for the second step
- Read the displayed temperature of the MEDIAWEL during the first step (after 5 min of stability at 121°C) and read the one of the reference thermometer at the same time and write them down
- Read the displayed temperature of the MEDIAWEL during the second step (after 5 min of stability at 42°C) and read the one of the reference thermometer at the same time and write them down
- Stop the cycle

THERMOMETER WITH PROBE  
ON REQUEST

- Access the service menu: “SERVICE”

The image shows two screenshots of the MEDIAWEL 10 interface. The left screenshot displays the password entry screen, featuring a lock icon, a 'Password' label, and a numeric keypad. The temperature display shows 19.2 °C. The right screenshot displays the calibration screen, with input fields for 'Real temperatures (reference probe)' and 'Measured temperatures (internal probe)', and a temperature display of 19.2 °C. A green checkmark button is visible in the bottom right corner of the right screenshot.

- Enter carefully in the top left box the temperature read on the reference thermometer during the second step
- Enter carefully in the top right box the temperature read on the reference thermometer during the first step
- Enter carefully in the bottom left box the temperature read on the reference thermometer during the second step
- Enter carefully in the bottom right box the temperature read on the reference thermometer during the first step.
- Make sure the data entered are correct and click on validate (green button)

The temperature on the main screen will adjust according to the calibration.