VIASURE

Real Time PCR Detection Kits

by **CerTest**
BIOTEC

Borrelia, Anaplasma & Coxiella

Handbook for the following references/

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VIASURE <i>Borrelia, Anaplasma & Coxiella</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-BAC106L
VIASURE <i>Borrelia, Anaplasma & Coxiella</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-BAC106H
VIASURE <i>Borrelia, Anaplasma & Coxiella</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-BAC112L
VIASURE <i>Borrelia, Anaplasma & Coxiella</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-BAC112H
VIASURE <i>Borrelia, Anaplasma & Coxiella</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-BAC113L
VIASURE <i>Borrelia, Anaplasma & Coxiella</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-BAC113H
VIASURE <i>Borrelia, Anaplasma & Coxiella</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-BAC136
VIASURE <i>Borrelia, Anaplasma & Coxiella</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-BAC172

ENGLISH

1. Intended use

VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Borrelia burgdorferi sensu lato (s.l.)*, *Borrelia miyamotoi* and/or *Borrelia hermsii*, *Anaplasma phagocitophylum* and/or *Coxiella burnetii* in clinical samples from patients with signs and symptoms of Tick-Borne diseases. This test is intended for use as an aid in the diagnosis of *B. burgdorferi sensu lato (s.l.)*, *B. miyamotoi* and/or *B. hermsii*, *Anaplasma phagocitophylum* and/or *Coxiella burnetii* in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, multiplied using Real Time amplification and detected using fluorescent reporter dye probes specific *Borrelia burgdorferi s.l.*, *Borrelia miyamotoi* and *Borrelia hermsii*, *Anaplasma phagocitophylum* and *Coxiella burnetii*.

2. Summary and Explanation

Tick Borne diseases comprise a group of infections transmitted to humans by the bite of ticks infected with bacteria, viruses, or parasites. Tick Borne diseases affecting humans include Lyme disease, relapsing fever, babesiosis, anaplasmosis, ehrlichiosis, Q fever, Tick Borne encephalitis and spotted fever among others.

Lyme disease (or Lyme borreliosis) is the most common Tick-Borne disease worldwide. It is caused by the spirochaete *Borrelia burgdorferi* and it is transmitted to humans through the bite of infected blacklegged ticks of the genus *Ixodes*. Typical symptoms of Lyme disease include fever, headache, fatigue and a characteristic skin rash called erythema migrans. If left untreated, infection can spread to joints, the heart, and the nervous system.

Borrelia miyamotoi and *B. hermsii* are species of spiral-shaped bacteria that is closely related to the bacteria that cause tick-borne relapsing fever (TBRF). First identified in 1995 in ticks from Japan, *B. miyamotoi* has also been detected in different *Ixodes* ticks species, whereas spirochete *Borrelia hermsii* is transmitted by its argasid tick vector, *Ornithodoros hermsi*. Patients with this infection were most likely to have fever, chills, and headache. Other common symptoms included body and joint pain and fatigue.

Anaplasmosis is caused by the bacterium *Anaplasma phagocytophilium* and is transmitted to humans through ticks of the genus *Ixodes*. Anaplasmosis shows high clinical variability, with symptoms like headache, fever, chills, malaise, muscle pain, nausea, cough, confusion and rash. If not treated correctly, anaplasmosis can evolve to severe clinical manifestations and even death in <1% of cases.

Q fever is a zoonosis caused by the bacteria *Coxiella burnetii*. Cattle, sheep, and goats are the primary reservoirs of *Coxiella burnetii*, and transmission to humans occurs primarily through inhalation of aerosols from contaminated soil or animal waste, though it can be transmitted through tick bites. Only about 50% of the infected people show clinical symptoms, which are flu-like symptoms such as headache, fever, chills, fatigue, muscle aches, nausea, cough, chest pain and weight loss. In severe cases people may develop pneumonia or hepatitis.

Since most Tick-Borne diseases show similar symptoms, diagnosis can be problematic. Real-time PCR assays have been shown to be a sensitive and specific diagnostic tool for the detection of the causative agent.

3. Principle of the procedure

VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit is designed for the diagnosis of *Borrelia burgdorferi s.l.*, *Borrelia miyamoto* and/or *Borrelia hermsii*, *Anaplasma phagocitophylum* and/or *Coxiella burnetii* in clinical samples. After DNA isolation, the identification of *Borrelia burgdorferi s.l.*, *Borrelia miyamoto* and/or *Borrelia hermsii*, *Anaplasma phagocitophylum* and/or *Coxiella burnetii* is performed by the amplification of a conserved region of the 23S rRNA gene (*Borrelia burgdorferi s.l.*, *Borrelia miyamotoi* and/or *Borrelia hermsii*), *msh2* gene (*Anaplasma phagocitophylum*) and *IS1111* gene (*Coxiella burnetii*), using specific primers and a fluorescent-labelled probe.

VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. *Borrelia burgdorferi s.l.*, *Borrelia miyamotoi* and *Borrelia hermsii* DNA targets are amplified and detected in FAM channel, *Anaplasma phagocitophylum* DNA targets are amplified and detected in ROX channel, and *Coxiella burnetii* DNA targets are amplified and detected in Cy5 channel. And the internal control (IC) which is amplified and detected in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>Borrelia, Anaplasma & Coxiella</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Borrelia, Anaplasma & Coxiella</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *Borrelia, Anaplasma & Coxiella* Real Time PCR Detection Kit with Ref. VS-BAC106L, VS-BAC106H, VS-BAC112L and VS-BAC112H.

Reagent/Material	Description	Colour	Amount
<i>Borrelia, Anaplasma & Coxiella</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Borrelia, Anaplasma & Coxiella</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *Borrelia, Anaplasma & Coxiella* Real Time PCR Detection Kit with Ref. VS-BAC113L, and VS-BAC113H.

Reagent/Material	Description	Colour	Amount
<i>Borrelia, Anaplasma & Coxiella</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Borrelia, Anaplasma & Coxiella</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 4-cap strip

Table 3. Reagents and materials provided in VIASURE *Borrelia, Anaplasma & Coxiella* Real Time PCR Detection Kit with Ref. VS-BAC136 and VS-BAC172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.

- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-BAC113L, VS-BAC113H, VS-BAC136 and VS-BAC172). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-BAC136 and VS-BAC172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from serum, blood and microbiological culture samples, you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions.

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- ZP02006 MagPurix Bacterial DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).

8.2. Lyophilized positive control

Borrelia, *Anaplasma* & *Coxiella* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the

lyophilized *Borrelia*, *Anaplasma* & *Coxiella* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA extracted from each sample, reconstituted *Borrelia*, *Anaplasma* & *Coxiella* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kits).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*Borrelia burgdorferi* s.l./*Borrelia miyamotoi*/*Borrelia hermsii*), HEX, JOE or VIC (Internal Control (IC)), ROX (*Anaplasma phagocitophilum*), and Cy5 (*Coxiella burnetii*) channels. Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005PT™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Borrelia*, *Anaplasma* & *Coxiella* in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

<i>Borrelia burgdorferi</i> s.l./ <i>Borrelia miyamotoi</i> / <i>Borrelia hermsii</i> (FAM)	<i>Anaplasma phagocitophyllum</i> (ROX)	<i>Coxiella burnetii</i> (Cy5)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	<i>B. burgdorferi</i> s.l./ <i>B. miyamotoi</i> / <i>B. hermsii</i> , <i>A. phagocitophyllum</i> and <i>C. burnetii</i> Positive
-	-	-	+	-	+	<i>B. burgdorferi</i> s.l./ <i>B. miyamotoi</i> / <i>B. hermsii</i> , <i>A. phagocitophyllum</i> and <i>C. burnetii</i> Negative
+	-	-	+/-	-	+	<i>B. burgdorferi</i> s.l./ <i>B. miyamotoi</i> / <i>B. hermsii</i> Positive, <i>A. phagocitophyllum</i> and <i>C. burnetii</i> Negative
+	+	-	+/-	-	+	<i>B. burgdorferi</i> s.l./ <i>B. miyamotoi</i> / <i>B. hermsii</i> and <i>A. phagocitophyllum</i> Positive, and <i>C. burnetii</i> Negative
+	-	+	+/-	-	+	<i>B. burgdorferi</i> s.l./ <i>B. miyamotoi</i> / <i>B. hermsii</i> and <i>C. burnetii</i> Positive, and <i>A. phagocitophyllum</i> Negative
-	+	-	+/-	-	+	<i>A. phagocitophyllum</i> Positive, <i>B. burgdorferi</i> s.l./ <i>B. miyamotoi</i> / <i>B. hermsii</i> and <i>C. burnetii</i> Negative
-	+	+	+/-	-	+	<i>A. phagocitophyllum</i> and <i>C. burnetii</i> Positive, <i>B. burgdorferi</i> s.l./ <i>B. miyamotoi</i> / <i>B. hermsii</i> Negative
-	-	+	+/-	-	+	<i>C. burnetii</i> Positive, <i>B. burgdorferi</i> s.l./ <i>B. miyamotoi</i> / <i>B. hermsii</i> and <i>A. phagocitophyllum</i> Negative
-	-	-	-	-	+	Experiment fail
+	+	+	+	+	-	Experiment fail

Table 5. Sample interpretation

+: Amplification curve


-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from blood, serum, tissue samples and and microbiological culture.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.

- There is a possibility of false positive results due to cross-contamination by contamination by *Borrelia*, *Anaplasma* & *Coxiella*, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit was tested using 95 DNA samples. A total of 17 well characterized *Borrelia* strains comprising 9 different *Borrelia burgdorferi* sensu lato genospecies, notably *B. japonica* (n=1), *B. burgdorferi* sensu stricto (n=2; B31 and PBr strains), *B. bavariensis* (n=1; PBi strain), *B. garinii* (n=5; PBr, PHei, PWudII, PRef and PLa strains), *B. bissettii* (n=1; PGeb strain), *B. afzelii* (n=2; PKo and PVPM strains), *B. lusitaniae* (n=1; Poti B2strain), *B. spielmanii* (n=1; PSigII strain), *B. valaisiana* (n=1; VS116 strain)) were included. Additionally, two relapsing fever control strains *B. hermsii* (n=1) and *B. miyamotoi* (n=1) were included as well as the potentially cross-reactive spirochaetes *Leishmania* spp. (n=2) and *Treponema* spp. (n=1) were tested. The VIASURE real-time PCR successfully detected all tested *Borrelia* sensu lato genospecies and the relapsing fever group strains *B. hermsii* and *B. miyamotoi*. No cross reactivity was observed with DNA from *Leptospira* and *Treponema* species for VIASURE assay.

VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit was evaluated with INSTAND *Coxiella burnetii* & *Bacillus anthracis* panels from 2017 and 17 additional tissue samples. The results were compared with the final EQA program reports or with those obtained by a commercial qPCR assay (EXOone *Coxiella burnetii* (EXOPOL)). All *Coxiella burnetii* positive samples from INSTAND program were detected and 15/17 clinical samples showed a positive result in the identification of *Coxiella burnetii*.

In conclusion, the results show a high sensitivity and specificity to detect *Borrelia*, *Anaplasma* and *Coxiella* using VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction (Figure 2, 3, 4).

Figure 2. Dilution series of *Borrelia burgdorferi*/*Borrelia miyamotoi*/ *B. hermsii* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM)


Figure 3. Dilution series of *Anaplasma phagocitophilum* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel ROX).



Figure 4. Dilution series of *Coxiella burnetii* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel Cy5).

12.3. Analytical specificity

The specificity of *Borrelia*, *Anaplasma* & *Coxiella* assay was confirmed by testing a panel consisting of different microorganisms representing the most common Vector Borne and tropical pathogens. No cross-reactivity was detected against any of the following microorganisms tested:

Cross-reactivity testing					
<i>Anaplasma marginale</i>	-	West Nile Virus (Heja)	-	Zika virus strain FB-GWUH-2016	-
<i>Bartonella henselae</i> strain Houston-1	-	West Nile Virus (NY99)	-	Zika Virus (African strain)	-
<i>Theileria annulata</i>	-	West Nile Virus (Ug37)	-	Zika Virus (Asian strain PF13/251013-18)	-
<i>Rickettsia conorii</i> (Morocann strain)	-	Yellow Fever Virus (17D strain)	-	Zika Virus strain 11474/16 (French Polynesia)	-
<i>Leptospira</i>	-	Yellow Fever Virus French Neurotropic	-	Zika Virus strain 11468/16(French Polynesia)	-
<i>Treponema phagedenis</i>	-	Chikungunya virus strain S27 Petersfield	-	Dengue 1 virus strain Hawaii	-
<i>Trypanosoma cruzi</i>	-	Chikungunya virus Martinique isolate	-	Dengue 2 virus strain New Guinea C	-
<i>Japanese encephalitis</i>	-	Chikungunya virus strain F24	-	Dengue 3 virus strain H87	-
<i>Plasmodium falciparum</i> (3D7 strain)	-	Chikungunya virus WHO IS (R91064)	-	Dengue 4 virus strain H241	-
Tick Borne Encephalitis Virus (TBEV)	-	St Louis Encephalitis virus	-		

Table 6. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *Borrelia*, *Anaplasma* & *Coxiella* Real Time PCR Detection Kit was evaluated against *Borrelia azfelic* (P-Ko/1984 and PVPM strains), *Borrelia bavariensis* (PBi strain), *Borrelia bissettii* (PGeb strain), *Borrelia bissettiae*, *Borrelia burgdorferi* strain IRS, *Borrelia burgdorferi sensu stricto* (B31 and PBre strains), *Borrelia garinii* (PHei, PWudll, PRef and PLa strains) and *Borrelia garinii* OspA Typ3 strains, *B. japonica*, *B. lusitanae* (Poti B2strain), *B. spielmanii* (PSigll strain), *Borrelia valaisiana* (VS116 strain), *B. hermsii* and *B. miyamotoi*, *Anaplasma phagocitophilum* and *Coxiella burnetii* strain Nine Mile Q showing positive results.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

(1) Select Ramp Speed "Standard".
 (2) See Annex 3 to check optical measurement exposure setting.
 (3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.
 (4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.
 (5) No detection in Cy5 channel.
 (6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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

13. Bibliography/Bibliografía

1. V. Parikh *et al.* Infections of the nervous system. *International Journal of Critical Illness and Injury Science* 2012; 2(2): 82–97.
2. J. Brites-Neto *et al.* Tick Borne infections in human and animal population worldwide. *Veterinary World* 2015, 8(3):301-315.
3. J.J. Sormunen *et al.* Tick-borne bacterial pathogens in southwestern Finland. *Parasites & Vectors*. 2016; 9: 168.
4. M.P. Nelder *et al.* Human pathogens associated with the blacklegged tick *Ixodes scapularis*: a systematic review. *Parasites & Vectors* 2016; 9:265
5. M.A. Diuk-Wasser *et al.* Coinfection by the Tick Borne pathogens *Babesia microti* and *Borrelia burgdorferi*: ecological, epidemiological and clinical consequences. *Trends in Parasitology* 2016; 32(1): 30–42.
6. A. Hojgaard *et al.* Detection of *Borrelia burgdorferi*, *Anaplasma phagocytophilum* and *Babesia microti*, with two different multiplex PCR assays. *Tick & Tick Borne Diseases* 2014; 5(3): 349-351.
7. J. Liu *et al.* Development of a TaqMan Array Card for Acute-Febrile-Illness, Outbreak Investigation and Surveillance of Emerging Pathogens, Including Ebola Virus. *Journal of Clinical Microbiology* 2016; 54(1): 49-58.

8. Centers for Disease Control and Prevention (<https://www.cdc.gov/niosh/topics/Tick Borne/>)
9. World Health Organization (<http://www.who.int/en/>).

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™, ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: February 2020

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev01

VIASURE

Real Time PCR Detection Kits

CerTest
BIOTEC