VIASURE

Real Time PCR Detection Kits

by CerTest
BIOTEC

FOR INFORMATION PURPOSES ONLY

Pneumocystis jirovecii

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VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-JIR106L
VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-JIR106H
VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-JIR112L
VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-JIR112H
VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-JIR113L
VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-JIR113H
VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-JIR136
VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-JIR172

ENGLISH

1. Intended use

VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit is designed for qualitative and quantitative detection of *Pneumocystis jirovecii* in respiratory samples from patients with signs and symptoms of respiratory infection. This test is intended for use as an aid in the diagnosis of *Pneumocystis jirovecii* in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *Pneumocystis jirovecii*.

2. Summary and Explanation

Pneumocystis jirovecii pneumonia (PCP) is an acute and life-threatening lung disease caused by the fungus *Pneumocystis jirovecii*. PCP is an important disease of immunocompromised humans, particularly patients with HIV, but also patients with an immune system that is severely suppressed for other reasons. In humans with a normal immune system, it is an extremely common silent infection. In developing regions of the world, the prevalence of PCP was once thought to be much lower, but studies have shown that the lower reported incidence is likely a failure to accurately diagnose.

The symptoms of PCP are nonspecific, in patients with HIV tends to present much later, often after several weeks of symptoms, compared with PCP associated with other immunocompromising conditions. Symptoms of PCP include the following: progressive exertional dyspnea, fever, non-productive cough, chest discomfort, weight loss, chills and hemoptysis (rare).

PCP is difficult to diagnose as a result of the associated nonspecific signs and symptoms. Because *P. jirovecii* cannot be propagated in culture, microscopic visualization of cysts or trophic forms in pulmonary specimens with cytochemical or immunofluorescent staining with monoclonal antibodies and/or DNA amplification are the standard procedures to detect this microorganism.

3. Principle of the procedure

VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit is designed for the diagnosis of *Pneumocystis jirovecii* in respiratory samples. After DNA isolation, the identification of *Pneumocystis jirovecii* is performed by the amplification of a conserved region of the large-subunit (mt LSU) rRNA gene using specific primers and a fluorescent-labelled probe.

VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on real time PCR platforms. In addition, quantification of specific *Pneumocystis jirovecii* DNA can be achieved by generating of a standard curve using the *Pneumocystis jirovecii* Quantitative Standard provided with the kit.

VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. *Pneumocystis jirovecii* DNA targets are amplified and detected in FAM channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>Pneumocystis jirovecii</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Pneumocystis jirovecii</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
<i>Pneumocystis jirovecii</i> Quantitative Standard	Quantitative Standard	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit with Ref. VS-JIR106L, VS-JIR106H, VS-JIR112L and VS-JIR112H.

Reagent/Material	Description	Color	Amount
<i>Pneumocystis jirovecii</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Pneumocystis jirovecii</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
<i>Pneumocystis jirovecii</i> Quantitative Standard	Quantitative Standard	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit with Ref VS-JIR113L and VS-JIR113H.

Reagent/Material	Description	Colour	Amount
<i>Pneumocystis jirovecii</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Pneumocystis jirovecii</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
<i>Pneumocystis jirovecii</i> Quantitative Standard	Quantitative Standard	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 4-cap strip

Table 3. Reagents and materials provided in VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit with Ref. VS-JIR136 and VS-JIR172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid) and Roche Molecular Diagnostics Cobas z480 Analyzer. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.

- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7 Precautions for users

- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-JIR113L, VS-JIR113H, VS-JIR136 and VS-JIR172). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-JIR136 and VS-JIR172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Do not confuse *Pneumocystis jirovecii* Positive Control vial with *Pneumocystis jirovecii* Quantitative Standard vial.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from respiratory samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- EZ1 Virus Mini Kit, using EZ1 instrument (Qiagen).

8.2. Lyophilized positive control and quantitative standard

Pneumocystis jirovecii Positive Control and *Pneumocystis jirovecii* Quantitative Standard contain high copies template, the recommendation is to open and manipulate them in a separate laboratory area away from the other components. Reconstitute the lyophilized *Pneumocystis jirovecii* Positive Control and/or *Pneumocystis jirovecii* Quantitative Standard (red vials) adding 100 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. Standard curve preparation for quantitative assay

This kit can be used to quantify the DNA copy number of *Pneumocystis jirovecii*. To perform a quantitative assay, a standard curve must be prepared by serial dilution of *Pneumocystis jirovecii* Quantitative Standard (which contains approximately 2×10^7 copies/µl*). Take *Pneumocystis jirovecii* Quantitative Standard (white pouch, red vial) as the starting high standard in the first tube, and prepare the standard curve dilution series as follows:


- Pipette 90 µl of RNase/DNase free Water in 6 microcentrifuge tubes (1.5 mL or 2 mL).
- Add 10 µl of *Pneumocystis jirovecii* Quantitative Standard to the first tube to get a standard with about 2×10^6 copies/µl. Mix by vortex and centrifuge briefly.
- Add 10 µl of standard with $\sim 2 \times 10^6$ copies/µl to the second tube to get a standard with approximately 2×10^5 copies/µl. Mix by vortex and centrifuge briefly.
- Repeat the previous step sequentially to complete the dilution series for standards with about 2×10^4 , 2×10^3 , 2×10^2 and 2×10^1 copies/µl (See Figure 1).

* For more details see "Certificate of Analysis", where lot number and DNA copy number of *Pneumocystis jirovecii* Quantitative Standard are detailed.

Dilution is not needed for performance of qualitative real-time PCR assay. In those assays, recommendation is to use *Pneumocystis jirovecii* Positive Control to minimize the risk of cross-contamination.

Figure 1. Preparation of the standard curve dilution series for quantitative assay

8.4. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, Negative Control (violet vial), *Pneumocystis jirovecii* Positive Control (aluminum pouch, red vial) for qualitative test or standard curve dilutions for quantitative assays in different wells, and close the wells with the caps provided.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kits).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*Pneumocystis jirovecii*) and HEX, JOE or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

9.1. Qualitative assay

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Pneumocystis jirovecii* in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

<i>Pneumocystis jirovecii</i> (FAM)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
+	+/-	-	+	<i>Pneumocystis jirovecii</i> Positive
-	+	-	+	<i>Pneumocystis jirovecii</i> Negative
-	-	-	-	Experiment fail
+	+	+	+	Experiment fail

Table 5. Sample interpretation

+: Amplification curve


-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 2. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

Negative control

Positive control

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

9.2. Quantitative assay

A standard curve can be generated from the dilution series of *Pneumocystis jirovecii* Quantitative Standard following the formula:

$$Ct = m \log (Q) + b$$

Where: Ct = Threshold Cycle

m = Slope

Q = Concentration


b = Intercept

Positive samples of unknown concentration can be quantified by interpolating their Ct value in the standard curve according to the formula:

$$Q = 10^{(Ct-b)/m}$$

Figure 3. Standard curve from the dilution series (blue) and quantification of a positive sample of unknown concentration (orange) run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The concentration of your "Positive sample" is displayed in copies/ μ l and refers to the concentration in the Eluted DNA not to the original clinical sample. To determine the target concentration of the original sample, take into account the dilutions of the extraction procedure and PCR Set-up

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from bronchoalveolar lavage (BAL) specimens.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Pneumocystis jirovecii*, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit was tested using 14 bronchoalveolar lavage from patients with suspicious of pneumocystis pneumonia. These results were compared with those obtained by other real-time PCR commercial kit (FTD *Pneumocystis jirovecii*, Fast-Track DIAGNOSTICS).

The results were as follows:

VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit	FTD <i>Pneumocystis jirovecii</i> (Fast-Track DIAGNOSTICS)		
		+	-
+	11	0	11
-	0	3	3
Total	11	3	14


Table 6. Comparative results.

The results show a high sensitivity and specificity to detect *Pneumocystis jirovecii* using VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction (Figure 4).

Figure 4. Dilution series of *Pneumocystis jirovecii* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

12.3. Analytical specificity

The specificity of the *Pneumocystis jirovecii* assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing					
<i>Staphylococcus aureus</i> subsp. <i>Aureus</i>	-	Methicillin-resistant <i>Staphylococcus aureus</i>	-	Influenza B/Brisbane/60/2008-like virus	-
<i>Bordetella pertussis</i>	-	Respiratory syncytial virus (RSV)	-	Influenza B/Florida/04/06 virus	-
<i>Fluoribacter bozemanæ</i>	-	<i>Haemophilus influenzae</i>	-	Influenza B/Phuket/3073/2013 virus	-
<i>Legionella micdadei</i>	-	<i>Moraxella catarrhalis</i>	-	Human Adenovirus	-
<i>Legionella dumoffii</i>	-	Influenza A/California/7/2009(H1N1) virus	-	Human coronavirus 229E	-
<i>Legionella longbeachæ</i>	-	Influenza A/Perth/16/2009(H3N2) virus	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
<i>Legionella pneumophila</i>	-	Influenza A/New Caledonia/20/99(H1N1) virus	-	Human metapneumovirus A and B	-
<i>Mycoplasma pneumoniae</i>	-	Influenza A/Switzerland/9715293/2013 virus	-	Human rhinovirus	-
<i>Streptococcus pneumoniae</i>	-	Influenza A/Turkey/Germany R2485+86/2014 virus	-		

Table 7. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit was evaluated against *Pneumocystis jirovecii* Type 1A, *Pneumocystis jirovecii* g885652 and *Pneumocystis jirovecii* j888023, showing positive results.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

(1) Select Ramp Speed "Standard".
 (2) See Annex 3 to check optical measurement exposure setting.
 (3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.
 (4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.
 (5) No detection in Cy5 channel.
 (6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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VIASURE <i>Pneumocystis jirovecii</i> Real Time PCR Detection Kit	FTD <i>Pneumocystis jirovecii</i> (Fast-Track DIAGNOSTICS)			
		+	-	Total
+		11	0	11
-		0	3	3
Total		11	3	14

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

13. Bibliography/Bibliografía

1. A. Roux *et al.* *Pneumocystis jirovecii* pneumonia in patients with or without AIDS, France. *Emerging Infectious Diseases* journal 2014; 20: 1490–1497.
2. E.J. Calderón *et al.* *Pneumocystis* infection in humans: diagnosis and treatment. *Expert Review of Anti-infective Therapy* 2010; 8: 683–701.
3. J.R. Harris *et al.* *Pneumocystis Jirovecii* Pneumonia: Current Knowledge and Outstanding Public Health Issues. *Current Fungal Infection Reports Journal*, 2010; 4(4): 229-237.
4. P. Rohner *et al.* Detection of *Pneumocystis jirovecii* by two staining methods and two quantitative PCR assays. *Infection*, 2009; 37(3):261-5.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: May 2019

FOR INFORMATION PURPOSES ONLY

FOR INFORMATION PURPOSES ONLY

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev01

VIASURE

Real Time PCR Detection Kits

CerTest
BIOTEC