# VIASURE

Real Time PCR Detection Kits

by CerTest  
BIOTEC

## *Campylobacter coli, C. lari & C. jejuni*

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---

### 1. Intended use

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Campylobacter coli*, *Campylobacter lari* and/or *Campylobacter jejuni* in human stool samples from patients with signs and symptoms of gastrointestinal infection. This test is intended to be used as an aid in the diagnosis of *C. coli*, *C. lari* and/or *C. jejuni* in combination with clinical and epidemiological risk factors. DNA is extracted from specimens, multiplied using Real Time amplification and detected using fluorescent reporter dye probes specific for *C. coli*, *C. lari* and *C. jejuni*.

### 2. Summary and Explanation

*Campylobacter* species are gram-negative, nonspore forming, spiral, or curved-shaped bacteria. Among the more than 26 species currently classified in the *Campylobacter* genus, most human diseases are attributed to three major food-borne species: *Campylobacter jejuni*, *Campylobacter coli*, and *Campylobacter lari*. These three species can be isolated from poultry and are of greatest concern to the poultry industry.

*Campylobacter* is considered one of the most common causes of diarrheal illness worldwide. Poultry is a major reservoir and source of transmission of *Campylobacter* to humans. In particular *C. jejuni* is the common species found in beef, *C. coli* is often isolated from pork, and *C. lari* is predominant in shorebirds. Other risk factors include consumption of animal products and water, contact with animals, and even person-to-person transmission (fecal-oral or via fomites).

Infection with *Campylobacter* causes gastroenteritis characterised by fever, vomiting, headaches, and abdominal pain with watery or bloody diarrhea, for a median duration of 6 days. Besides gastroenteritis, these three species can cause periodontitis, septicemia, and second trimester intrauterine growth restriction. Furthermore, *C. jejuni* infection may lead to autoimmune conditions such as Guillain-Barré syndrome (GBS) and Miller Fisher syndrome (MFS).

*Campylobacter* enteritis is usually self-limiting and typically does not require antimicrobial therapy. In these cases, maintenance of proper hydration and electrolyte balance is the most important tenets of treatment. However, in severe and prolonged cases of enteritis, bacteremia, or other extraintestinal infection, prompt antimicrobial treatment is indicated.

Traditional microbiological methods for *Campylobacter* identification include enrichment, culturing, isolation, and phenotypic characterization. The procedures are laborintensive, time consuming, and with a relatively narrow differentiation spectrum among target species. These factors present challenges for the identification of *Campylobacter* from patient samples or contaminated food. Fortunately, methodologies based on molecular biology have been developed to improve laboratory approaches, such as Real Time PCR. Multiplex qPCR can detect several targets in one reaction, saving time, effort, and sample. This qPCR assay can be used for identifying *C. coli*, *C. lari* and *C. jejuni* or secondary screening for confirmation of *Campylobacter* bacteria to the species level.

### 3. Principle of the procedure

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit is designed for the diagnosis of *H. influenzae*, *N. meningitidis* and/or *S. pneumoniae* in clinical samples. After DNA isolation, the identification of *C. coli*, *C. lari* and *C. jejuni* is performed by the amplification of a conserved region of the *hipO* gene for *Campylobacter jejuni*, *Gyrasa A* gene for *Campylobacter lari* and *CeuE* gene for *Campylobacter coli*, using specific primers and a fluorescent-labelled probe.

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format. The assay can use an Extraction Control (EC) which can be introduced into each sample at the lysis buffer stage of the extraction process. This control can be used to monitor the extraction process and/or discard the inhibition of the polymerase activity. *Campylobacter jejuni* DNA targets are amplified and detected in the FAM channel, *Campylobacter lari* DNA targets are amplified and detected in the ROX channel, *Campylobacter coli* DNA targets are amplified and detected in the Cy5 channel and the extraction control (EC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

### 4. Reagents provided

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table 1 and Table 2:

Reference	Reagent/Material	Description	Colour	Amount
VS-1CLJ1SL/ VS-1CLJ1SH	<i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	6/12 x 8-well strip
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-CLJ1C	<i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
VS-EC01	Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit with Ref. VS-CLJ106LE, VS-CLJ106HE, VS-CLJ112LE and VS-CLJ112HE.

Reference	Reagent/Material	Description	Color	Amount
VS-1CLJ1PL/ VS-1CLJ1PH	<i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-CLJ1C	<i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
VS-EC01	Extraction Control	Non-infectious nucleic acid lyophilized	Green	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit with Ref VS-CLJ113LE and VS-CLJ113HE.

## 5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.

- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

## 7. Precautions for users

- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-CLJ113LE and VS-CLJ113HE). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.

## 8. Test procedure

### 8.1. Lyophilized extraction control

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the lyophilized positive control. Reconstitute the lyophilized Extraction Control (green vial) adding 500 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C. Recommendation is to separate it in aliquots to minimize freeze and thaw cycles.

Note: The Water RNase/DNase free vial must be utilized first to reconstitute the lyophilized Extraction Control in pre-PCR laboratory area, and subsequently, it can be used for reconstitute the lyophilized *Campylobacter coli*, *C. lari* & *C. jejuni* Positive Control in an area away from the other components.

## 8.2. Sample preparation

Stool samples should be collected in clean containers and processed as soon as possible to guarantee the quality of the test. We recommend to use fresh samples.

For longer storage, the samples must be frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation. Freezing and thawing cycles are not recommended.

Perform the sample preparation following the manual system protocol or according to the recommendations appearing in the instructions for use of extraction kit used.

### 8.2.1. DNA extraction

If the Extraction Control is used to monitor nucleic acid isolation and as PCR inhibition control, add 5 µl of the EC to the specimen-lysis buffer mixture. The EC must not be added directly to the sample. Close each tube and vortex for 10 seconds.

If the Extraction Control is used only as a PCR inhibition control, 1µl of the EC should be added to the reconstituted Reaction-Mix.

For DNA extraction from stool samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Invisorb® Spin Universal Kit (Strattec).
- NucleoMag® Pathogen (Macherey Nagel).
- NucleoSpin® RNA Virus (Macherey Nagel).
- Maxwell® RSC Blood DNA Kit, using the Maxwell® 16 instrument (Promega).
- MagDEA Dx SV kit, using the magLEAD® 6gC instrument (Precision System Science Co.)
- NX-48 Stool DNA Kit, using the Nextractor® NX-48 system (Genolution).
- ZP02011 MagPurix Viral/Pathogen Nucleic Acids Extraction Kit A, using the MagPurix 12A instrument (Zinexts Life Science Corp.).
- ZP02012 MagPurix Viral/Pathogen Nucleic Acids Extraction Kit B, using the MagPurix 12A instrument (Zinexts Life Science Corp.).

## 8.3. Lyophilized positive control

*Campylobacter coli*, *C. lari* & *C. jejuni* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the

lyophilized *Campylobacter coli*, *C. lari* & *C. jejuni* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

### 8.4. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of Negative Control (violet vial) and 1 µL of the EC (green vial) in the reserved wells for negative control. The EC is used as extraction and/or PCR inhibition control.

Add 5 µL of DNA extracted from each sample in different wells. If the EC is only used as PCR inhibition control, we recommend to add 1 µL of the EC (green vial) to the sample wells.

Add 5 µL of reconstituted *Campylobacter coli*, *C. lari* & *C. jejuni* Positive Control (red vial) and 1 µL of the EC (green vial) in the reserved wells for positive control. The EC is used as extraction and/or PCR inhibition control.

Close the wells with the caps provided. It is recommended to briefly centrifuge the 8-well strips or 96-well plate.

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 3. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the FAM (*C. jejuni*), ROX (*C. lari*), Cy5 (*C. coli*) and HEX, JOE or VIC channels (Extraction Control (EC)). Depending on the equipment used select the proper detection channel (see Annex 2). In the Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

## 9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Campylobacter coli*, *C. lari* & *C. jejuni* in the positive control well. Check Extraction Control (EC) signal to verify the extraction procedure and correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

<i>C. jejuni</i> (FAM)	<i>C. lari</i> (ROX)	<i>C. coli</i> (Cy5)	Extraction control (HEX)	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	<i>C. jejuni</i> , <i>C. lari</i> and <i>C. coli</i> Positive
-	-	-	+	-	+	<i>C. jejuni</i> , <i>C. lari</i> and <i>C. coli</i> Negative
+	-	-	+/-	-	+	<i>C. jejuni</i> Positive, <i>C. lari</i> and <i>C. coli</i> Negative
+	+	-	+/-	-	+	<i>C. jejuni</i> and <i>C. lari</i> Positive, and <i>C. coli</i> Negative
+	-	+	+/-	-	+	<i>C. jejuni</i> and <i>C. coli</i> Positive, and <i>C. lari</i> Negative
-	+	-	+/-	-	+	<i>C. lari</i> Positive, <i>C. jejuni</i> and <i>C. coli</i> Negative
-	+	+	+/-	-	+	<i>C. lari</i> and <i>C. coli</i> Positive, <i>C. jejuni</i> Negative
-	-	+	+/-	-	+	<i>C. coli</i> Positive, <i>C. jejuni</i> and <i>C. lari</i> Negative
-	-	-	-	-	+	Experiment fail
+	+	+	+	+	-	Experiment fail

Table 4. Sample interpretation

+: Amplification curve

-: No amplification curve


A sample is considered positive if the Ct value obtained is less than 40 and the Extraction Control shows or not an amplification signal. Sometimes, the detection of Extraction Control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the Extraction Control is positive. A failure in the extraction procedure and/or an inhibition of the PCR reaction can be excluded by the amplification of Extraction Control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of signal in both, extraction control and target in sample wells, we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of nucleic acids purification and/or inhibition.

## 10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from stool samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Campylobacter coli*, *Campylobacter lari* and *Campylobacter jejuni* either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- There is a possibility of false positive results due to cross-contamination between Extraction Control and *Campylobacter coli*, *Campylobacter lari* and *Campylobacter jejuni* Positive Control, which contains high copies template, during their reconstitution by adding of Water RNase/DNase free (white vial). Each procedure must take place in established order and in a separate laboratory areas.

## 11. Quality control

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the Extraction Control (EC) in each well confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. Clinical sensitivity and specificity

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit was tested using 62 faecal samples (for *C. lari* are spiked clinical samples) from symptomatic patients. These results were compared with those obtained by Routine Diagnostic procedures (Microbiological culture and MALDI-TOF).

The results were as follows:

VIASURE <i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Real Time PCR Detection Kit	Culture + MALDI-TOF		
		+	-
+	10	0	10
-	2	50	52
Total	12	50	62

Table 5. Comparative results for *C. coli*

VIASURE <i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Real Time PCR Detection Kit	Culture + MALDI-TOF		
		+	-
+	30	0	30
-	3	29	32
Total	33	29	62

Table 6. Comparative results for *C. jejuni*

VIASURE <i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Real Time PCR Detection Kit	Culture + MALDI-TOF		
		+	-
+	5	0	5
-	0	57	57
Total	5	57	62

Table 7. Comparative results for *C. lari*

The results show a high sensitivity and specificity to detect *C. coli*, *C. lari* and *C. jejuni* using VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit.

## 12.2. Analytical sensitivity

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit has a detection limit of  $\geq 10$  DNA copies per reaction for *Campylobacter coli*, *Campylobacter lari* and *Campylobacter jejuni* (Figure 2, 3 and 4).

Figure 2. Dilution series of *Campylobacter jejuni* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).


Figure 3. Dilution series of *Campylobacter lari* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel ROX).



Figure 4. Dilution series of *Campylobacter coli* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel Cy5).

### 12.3. Analytical specificity

The specificity of the *Campylobacter coli*, *Campylobacter lari* and *Campylobacter jejuni* assay was confirmed by testing a panel consisting of different microorganisms representing the most common enteric pathogens or flora present in the intestine. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>Shigella flexneri</i>	-	<i>Campylobacter hyointestinalis</i>	-	<i>Dientamoba fragilis</i>	-
<i>Shigella dysenteriae</i> serotype 1	-	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	-/+	<i>Enterococcus faecalis</i>	-
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>typhi</i>	-	<i>Campylobacter lari</i>	-/+	<i>Blastocystis hominis</i>	-
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>paratyphi A</i>	-	<i>Campylobacter upsaliensis</i>	-	<i>Giardia intestinalis</i>	-
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>paratyphi B</i>	-	<i>Citrobacter freundii</i>	-	<i>Yersinia enterocolitica</i> O:3	-
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>typhimurium</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	<i>Yersinia enterocolitica</i> O:9	-
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>bongori</i> serovar 66:z41	-	<i>Serratia liquefaciens</i>	-	<i>Bacteroides fragilis</i>	-
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>enteritidis</i>	-	<i>Vibrio parahaemolyticus</i> serotype O1:K1	-	<i>Cryptosporidium parvum/hominis</i>	-
<i>Salmonella pullorum</i>	-	<i>Clostridium difficile</i>	-	<i>Entamoeba dispar</i>	-
<i>Salmonella gallinarum</i>	-	<i>Clostridium perfringens</i>	-	<i>Entamoeba histolytica</i>	-
<i>Helicobacter pylori</i> J99	-	<i>Escherichia coli</i> enterohemorrhagic serotype O157:H7	-	Sapovirus	-
<i>Helicobacter hepaticus</i>	-	<i>Escherichia coli</i> enteroinvasive	-	Rotavirus A	-
<i>Helicobacter cinaedi</i>	-	Enterotoxigenic <i>Escherichia coli</i> serotype O25:H42	-	Norovirus Genotypes I and II	-
<i>Helicobacter heilmannii</i>	-	Enteropathogenic <i>Escherichia coli</i>	-	Astrovirus Genotype I-VIII	-
<i>Proteus vulgaris</i>	-	<i>Klebsiella oxytoca</i>	-	Adenovirus types I-V	-
<i>Aeromonas caviae</i>	-	<i>Listeria monocytogenes</i> serovar 1/2c	-	Adenovirus type VIII and XV	-
<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i>	-	<i>Candida albicans</i>	-	Adenovirus serotype 40	-
<i>Campylobacter coli</i>	-/+	<i>Arcobacter butzleri</i>	-	Adenovirus serotype 41	-
<i>Campylobacter fetus</i> subsp. <i>fetus</i>	-	<i>Pseudomonas aeruginosa</i>	-		

Table 7. Reference pathogenic microorganisms used in this study.

### 12.4. Analytical reactivity

The reactivity of VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit for *Campylobacter coli* was evaluated against *Campylobacter coli* showing positive result.

The reactivity of VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit for *Campylobacter lari* was evaluated against *Campylobacter lari* showing positive result.

The reactivity of VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit for *Campylobacter jejuni* was evaluated against *Campylobacter jejuni* showing positive result.

ANNEX 1

**COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT**

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	7500 Fast Dx Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	StepOne™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System <sup>(6)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Roche	LightCycler® 480 Real-Time PCR System <sup>(4)</sup>
Roche	LightCycler® 96 Real-Time PCR System <sup>(4)</sup>
Roche	Cobas z480 Analyzer <sup>(4)</sup>

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ABI PRISM 7000 <sup>(6)</sup>
Applied Biosystems	ABI PRISM 7700 <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System <sup>(6)</sup>
Bio-Rad	MyiQ™2 Real-Time PCR Detection System <sup>(6)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
DNA-Technology	DTprime Real-time Detection Thermal Cycler <sup>(2)</sup>
DNA-Technology	DTlite Real-Time PCR System <sup>(2)</sup>
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System <sup>(2)</sup>
VIASURE	VIASURE 96 Real Time PCR System <sup>(2)</sup>

(1) Select Ramp Speed "Standard".  
 (2) See Annex 3 to check optical measurement exposure setting.  
 (3) The product should be reconstituted following the appropriate procedure (see Test Procedure, section 8.3) and transferred into the specific tubes designed to perform on Rotor-Gene® Q or SmartCycler® instruments.  
 (4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.  
 (5) No detection in Cy5 channel.  
 (6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

## ANNEX 2

## DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

## ANNEX 3

**OPTICAL MEASUREMENT EXPOSURE SETTING**

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits".

Universal exposition values are as follow:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -1000, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel – 500.

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
### 13. Bibliography/Bibliografía

1. N.O. Kaakoush *et al.* Global Epidemiology of *Campylobacter* Infection. *Clinical Microbiology Reviews* 2015; 28(3): 687-720.
2. S.M. Man. The clinical importance of emerging *Campylobacter* species. *Nature Reviews Gastroenterology & Hepatology* 2011; 8(12): 669-685.
3. R.F. de Boer *et al.* Detection of *Campylobacter* species and *Arcobacter butzleri* in stool samples by use of real-time multiplex PCR. *Journal of Clinical Microbiology* 2013; 51(1): 253-259.
4. K. C. Liu *et al.* Simultaneous Identification of *Campylobacter jejuni*, *Campylobacter coli*, and *Campylobacter lari* with SmartCycler-Based Multiplex Quantitative Polymerase Chain Reaction. *Foodborne pathogens and disease*, 2017 Volume 14, Number 7.
5. M. E. Patrick *et al.* Features of illnesses caused by five species of *Campylobacter*, Foodborne Diseases Active Surveillance Network (FoodNet) – 2010–2015. *Epidemiology & Infection* 2017. doi:10.1017/S0950268817002370.
6. C. A. Whitehouse *et al.* Antimicrobial Resistance in *Campylobacter* Species: Mechanisms and Genomic Epidemiology. *Advances in Applied Microbiology*, 2018, Volume 103.
7. AM. Mayr *et al.* Rapid detection and differentiation of *Campylobacter jejuni*, *Campylobacter coli*, and *Campylobacter lari* in food, using multiplex real-time PCR. *Journal of Food Protection*, 2010.73(2):241-50.

### 14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: August 2018

**CerTest Biotec, S.L.**

Pol. Industrial Río Gállego II - Calle J, Nº1  
50840, San Mateo de Gállego, Zaragoza (Spain)  
[www.certest.es](http://www.certest.es)

VIASURE online

F-362 rev01

**VIASURE**

Real Time PCR Detection Kits

**CerTest**  
BIOTEC