

RANDOX

INFECTIOUS DISEASE CONTROLS
(SEROLOGY)



ACUSERA



INFECTIOUS DISEASE CONTROLS

Third party, multi-marker controls designed to deliver accurate & reliable results



01 BENEFITS

02 INFECTIOUS DISEASE CONTROLS

03 CONTACT US

BENEFITS

Our Acusera Infectious Disease (Serology) Controls are designed to deliver a cost effective, high quality solution for the analysis of infectious diseases whilst producing trustworthy results time and time again. Covering a wide range of infectious diseases, including HIV, Hepatitis, EBV, ToRCH, and Lyme Disease, our controls are compatible for use on some of the most popular immunoassay platforms including Roche, Abbott, Siemens and Beckman among others.



Liquid ready-to-use

- All samples are provided in a user-friendly, liquid ready-to-use format significantly reducing preparation time and the risk of pipetting errors.



Commutability

- Manufactured from human plasma each control is designed to react to the test system in the same manner as the patient sample, helping to meet ISO 15189:2012 requirements.



True third party controls

- Designed to deliver an independent, unbiased assessment of performance with any instrument or method the Acusera range will help to meet ISO 15189:2012 requirements and eliminate the need for multiple instrument specific controls.



Reduced waste

- With a working stability of 60 days at 2°C to 8°C waste and costs are kept to a minimum.



Consolidation

- The availability of multi-marker controls will reduce the number of individual controls required ultimately reducing costs, time and storage space.



Clinically relevant

- Designed to challenge clinically relevant ranges our low reactive samples will ensure accurate and reliable instrument performance.



Acusera 24•7

- Our Acusera Infectious Disease (Serology) Controls are compatible for use with our interlaboratory data management software Acusera 24•7.

Lyme Disease (*Borrelia burgdorferi*)

Lyme Disease is an infectious disease caused by the bacterium *Borrelia burgdorferi* and spread by ticks.

Our control delivers a true third-party solution for the detection of Lyme Disease on most immunoassay analysers. All samples are conveniently supplied in a user-friendly, liquid ready-to-use format.

Analytes

Borrelia burgdorferi IgG | *Borrelia burgdorferi* IgM

Description

Description	Size	Cat. No.
Lyme Disease Negative Control	1 x 5 ml	SR10345
Lyme Disease Positive Control	1 x 5 ml	SR10346

Benefits

- Liquid ready-to-use
- Human based plasma
- Open vial stability of 60 days at 2°C to 8°C
- Suitable for use with most immunoassay analysers

ToRCH

ToRCH is most commonly associated with congenital anomalies caused by a variety of pathogens including *Treponema pallidum* (syphilis), Varicella Zoster Virus, Rubella, and Cytomegalovirus. Most ToRCH infections cause mild maternal morbidity but can have serious consequences on the foetus.

Our ToRCH portfolio includes positive controls for both IgG and IgM antibodies in addition to a negative control. Each control is manufactured using human plasma and is suitable for use with most immunoassay analysers. The availability of liquid ready-to-use samples helps to reduce preparation time and the potential for human error.

Analytes - ToRCH Negative

Cytomegalovirus (CMV) IgG | Cytomegalovirus (CMV) IgM | Epstein Barr Virus (EBV) EBNA IgG | Epstein Barr Virus (EBV) VCA IgG
 Epstein Barr Virus (EBV) IgM | *Helicobacter pylori* IgG | Herpes Simplex Virus 1 (HSV-1) IgG | Herpes Simplex Virus 1 (HSV-1) IgM
 Herpes Simplex Virus 2 (HSV-2) IgG | Herpes Simplex Virus 2 (HSV-2) IgM | Measles IgG | Mumps IgG | Rubella IgG | Rubella IgM
 Toxoplasma *gondii* IgG | Toxoplasma *gondii* IgM | *Treponema pallidum* (Syphilis) IgG | Varicella Zoster Virus (VZV) IgG

Analytes - ToRCH IgG Positive

Cytomegalovirus (CMV) IgG | *Helicobacter pylori* IgG | Herpes Simplex Virus 1 (HSV-1) IgG | Herpes Simplex Virus 2 (HSV-2) IgG
 Measles IgG | Mumps IgG | Rubella IgG | Toxoplasma *gondii* IgG | *Treponema pallidum* (Syphilis) IgG | Varicella Zoster Virus (VZV) IgG

Analytes - ToRCH IgM Positive

Cytomegalovirus (CMV) IgM | Herpes Simplex Virus 1 (HSV-1) IgM | Herpes Simplex Virus 2 (HSV-2) IgM | Rubella IgM
 Toxoplasma *gondii* IgM

Description

Description	Size	Cat. No.
ToRCH Negative Control	6 x 5 ml	SR10347
ToRCH IgG Positive Control	3 x 5 ml	SR10348
ToRCH IgM Positive Control	3 x 5 ml	SR10349

Benefits

- Liquid ready-to-use
- Human based serum
- Open vial stability of 60 days at 2°C to 8°C
- Suitable for use with most immunoassay analysers

INFECTIOUS DISEASE CONTROLS

Epstein Barr Virus (EBV)

EBV is an extremely common virus, spread primarily through saliva causing mononucleosis.

The Acusera EBV control is conveniently supplied liquid ready-to-use and is suitable for use with most immunoassay analysers.

Analytes

Epstein Barr Virus (EBV) EBNA IgG | Epstein Barr Virus (EBV) VCA IgG | Epstein Barr Virus (EBV) IgM

Description

EBV Positive Control

Size

1 x 5 ml

Cat. No.

SR10350

Benefits

- Liquid ready-to-use
- Human based serum
- Open vial stability of 60 days at 2°C to 8°C
- Suitable for use with most immunoassay analysers

Serology

A range of positive and negative serology controls are available comprising various infectious diseases including HIV, Hepatitis, HTLV and Syphilis.

Acusera Serology controls are supplied liquid ready-to-use and are suitable for use on most immunoassay analysers.

Analytes - Serology Negative

Anti-HAV | Anti-HBc | Anti-Hbe | Anti-HBs | Anti-HCV | Anti-HIV I / 2 | Anti-HTLV I / 2 | HAV IgM | HBc IgM | HbeAg | HBsAg
HIV P24Ag | *Treponema pallidum* (Syphilis) IgG

Analytes - Serology I Positive

Anti-HBc | Anti-HCV | Anti-HIV I / 2 | Anti-HTLV I / 2 | HBsAg | *Treponema pallidum* (Syphilis) IgG

Analytes - Serology II Positive

Anti-HAV | Anti-HBc | Anti-Hbe | Anti-HBs

Analytes - Serology III Positive

HAV IgM | HBc IgM

Description

Serology Negative Control

Size

6 x 5 ml

Cat. No.

SR10351

Serology I Positive Control

3 x 5 ml

SR10352

Serology II Positive Control

3 x 5 ml

SR10353

Serology III Positive Control

3 x 5 ml

SR10354

Benefits

- Liquid ready-to-use
- Human based plasma
- Open vial stability of 60 days at 2°C to 8°C
- Suitable for use with most immunoassay analysers

Contact us for more information on any of our products and services:

Headquarters

Randox Laboratories Ltd, 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom

☎ +44 (0) 28 9445 4399 ✉ marketing@randox.com 🌐 riqas.com

International Offices



AUSTRALIA

Randox (Australia) Pty Ltd.
Tel: +61 (0) 2 9615 4640



BRAZIL

Randox Brasil Ltda.
Tel: +55 11 5181-2024



CHINA

Randox Laboratories Ltd.
Tel: +86 021 6288 6240



CZECH REPUBLIC

Randox Laboratories S.R.O.
Tel: +420 2 1115 1661



FRANCE

Laboratoires Randox
Tel: +33 (0) 130 18 96 80



GERMANY

Randox Laboratories GmbH
Tel: +49 (0) 215 1937 0611



HONG KONG

Randox Laboratories Hong Kong Limited
Tel: +852 3595 0515



ITALY

Randox Laboratories Ltd.
Tel: +39 06 9896 8954



INDIA

Randox Laboratories India Pvt Ltd.
Tel: +91 80 6751 5000



POLAND

Randox Laboratories Polska Sp. z o.o.
Tel: +48 22 862 1080



PORTUGAL

Irandox Laboratorios Quimica Analitica Ltda
Tel: +351 22 589 8320



PUERTO RICO

Clinical Diagnostics of Puerto Rico, LLC
Tel: +1 787 701 7000



REPUBLIC OF IRELAND

Randox Teoranta
Tel: +353 7495 22600



SLOVAKIA

Randox S.R.O.
Tel: +421 2 6381 3324



SOUTH AFRICA

Randox Laboratories SA (Pty) Ltd.
Tel: +27 (0) 11 312 3590



SOUTH KOREA

Randox Korea
Tel: +82 (0) 31 478 3121



SPAIN

Laboratorios Randox S.L.
Tel: +34 93 475 09 64



SWITZERLAND

Randox Laboratories Ltd. (Switzerland)
Tel: +41 41 810 48 89



UAE

Randox Medical Equipments Trading LLC
Tel: +971 55 474 9075



USA

Randox Laboratories-US, Ltd.
Tel: +1 304 728 2890



VIETNAM

Randox Laboratories Ltd. Vietnam
Tel: +84 (0) 8 3911 0904

For technical support contact:

technical.services@randox.com




☎ +43 (0)2244 296 31

✉ info@labraum.at

🌐 www.labraum.at

RANDOX
QUALITY CONTROL



LT708 AUG20